

## **Mullard halogen quenched Geiger-Müller tube guide / Mullard Ltd.**

### **Contributors**

Mullard Limited. X-Ray Division.

### **Publication/Creation**

Mitcham Junction : Mullard, [ca.1960]

### **Persistent URL**

<https://wellcomecollection.org/works/hmu8fhpa>

### **License and attribution**

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome  
collection**

Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

**Mullard**

**halogen**

**quenched**

**Geiger-**

**Müller**

**tube**

**Guide**



# MULLARD GEIGER-MÜLLER TUBES

Mullard halogen quenched Geiger-Müller tubes have been designed to give years of satisfactory service under even the most exacting conditions. The high degree of mechanical and electrical strength, the exceptional stability, and the unique standardisation of characteristics have been achieved after a long period of intensive research into the selection of the most suitable materials and their processing.

Wherever possible the chrome line outside has been used as the main body of the tube, and the amount of glass required for insulation purposes has been reduced to a minimum. Most windows are sealed to the outside with a thin layer of glass, thus avoiding synthetic materials and the hazards associated with their use.

With the exception of the argon and krypton filled tubes for the measurement of X-radiation, all Mullard Geiger-Müller tubes are filled with neon, argon and nitrogen. A special technique has been developed for this process ensuring an exact and consistent specification which does not vary from tube to tube. The standardisation of working voltage is such that tubes taken at random or even supplied at different times will all operate under identical conditions. The minimum count rate is 10<sup>4</sup> and the shelf life is in excess of two years. Operation and life expectancy are not impaired by a moderate overvoltage, or by applying an operating voltage of incorrect polarity.

Type of counter	GAMMA			BETA END WINDOW		ALPHA BETA END WINDOW			X-RAY END WINDOW		HIGH CURRENT	LIQUID SAMPLE		Type of counter
Type number	MX112	MX115**	MX118	MX130**	MX113	MX122	MX118	MX122	MX119**	MX124	MX124 01			Type number
Threshold voltage (max)	300V	300V	300V	300V	300V	600V	1,000V	900V	900V	400V	400V			Threshold voltage (max)
Fluorescence length (max)	100V	100V	80V	100V	100V	100V	800V	200V	300V	100V	100V			Fluorescence length (max)
Fluorescence slope (average)	1%, 100V	0%, 100V	2%, 100V	1%, 100V	6%, 100V	1%, 100V	0%, 100V	0%, 100V	1%, 100V	1%, 100V	6%, 100V			Fluorescence slope (average)
Temp. range	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.	-65° C. to 75° C.			Temp. range
Dead time	140 µ sec.	100 µ sec.	110 µ sec.	110 µ sec.	80 µ sec.	80 µ sec.	100 µ sec.	100 µ sec.	35 µ sec.	100 µ sec.	100 µ sec.			Dead time
*Background (excluded)	30c/min.	20c/min.	50c/min.	21c/min.	6c/min.	70c/min.	90c/min.	50c/min.	30c/min.	30c/min.	30c/min.			*Background (excluded)
Background (included)	160c/min.	45c/min.	90c/min.	45c/min.	3.5 to 4.5 mg/cm <sup>2</sup>	1.6 to 2.5 mg/cm <sup>2</sup>	3.8 to 4.8 mg/cm <sup>2</sup>	3.5 to 4.5 mg/cm <sup>2</sup>						Background (included)
Window thickness														Window thickness
Window area														Window area
Wall thickness	400 mg/cm <sup>2</sup> *	315 mg/cm <sup>2</sup> *	375 mg/cm <sup>2</sup> *	222 mm <sup>2</sup>	31.7 mm <sup>2</sup> *	450 mm <sup>2</sup> *	314 mm <sup>2</sup> *	704 mm <sup>2</sup> *						Wall thickness
Active length	200 mm.	38 mm.	140 mm.											Active length
Minimum count rate	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	10 <sup>4</sup>	1,500 counts at 90 µA.	10 <sup>4</sup>	10 <sup>4</sup>			Minimum count rate
Electrical connections	Anode turret lug	Special base	Anode turret lug	Special base	Anode cap	Anode turret lug	Anode Cap	Anode Cap	Tinned leads	For use in mercury pools	For use in mercury pools			Electrical connections
Remarks	An all metal gas-tube for use in every measurement where high accuracy and identical mechanical robustness are required. Such applications include toxic, acidic, gaseous, and cosmic radiation measurements. Having rotator diameter, they may be used together in multiple arrangements.	An all metal tube in a separate socket, widely used in portable instruments. The tube may therefore be mounted externally on the equipment. The special socket is available separately if required—SO136.	A general purpose all metal gas-tube, which, owing to its greater length, is three times more sensitive than type MX115.	A low voltage metal end window tube which is interchangeable with tube MX115 when tube measurements are required. A metal window guard is provided under 100% nitrogen shield, the window area is 114 cm <sup>2</sup> but it may be removed if the full window area has to be used. Socket type MX130 is also used with this tube.	A miniature metal end window tube for all applications in which a very thin window of metal is required. These include in vivo measurements in the medical field, and scanning techniques where full advantage of the directional properties may be taken.	This versatile tube is carefully manufactured to exceptionally close tolerances in order that it may fit the majority of lead carrier and holders. The ultra thin metal window and low background make it suitable for the measurement of beta sources down to very low specific activities and energies.	A metal end window, argon-filled tube with a path length of approximately 10 centimetres for the measurement of X-radiation of wavelengths between 1.8 and 2.5 Angstroms. In fluorescence analysis it is suitable for elements of mass numbers 19 to 34.	A krypton-filled version of type MX118 for maximum efficiency in the measurement of X-radiation intensity by a micro-anne. It is suitable for the measurement of gamma radiation, or hard X-rays.	The tube passes a particularly high current (90 µA) at a radiation field of 10 r/hour for the direct reduction of radiation intensity by a micro-anne. It is suitable for the measurement of gamma radiation, or hard X-rays.	A liquid sample tube of outstanding stability throughout its life. Its liquid capacity is 0.10 ml. and a ground glass stopper (B01) is fitted to facilitate the application of solutions and the measurement of volatile fluids.	Type number MX124 01 is similar in all respects to MX124, but is not supplied with a stopper.			Remarks

\* Shielded by 10 mm. of lead plus 0.5 mm. aluminium.

\*\* Developed in collaboration with A.E.R.E. Harwell.

HALOGEN QUENCHED · ROBUST · LONG LIFE · RELIABLE · STANDARDISED CHARACTERISTICS

## Reliability . . .

. . . largely depends on the careful selection and use of components. This folder is a guide to the most reliable Geiger-Müller tubes available.

Mullard halogen quenched tubes are in a class apart for strength, dependability and consistent quality. There is never any need to order matched sets, as tubes of any one type operate at the same working voltage.

They withstand reasonable over-voltage, incorrect polarity, and, in the case of all metal types, surprisingly severe vibration and drop tests.

They are very long lived, and may be stored for years without deterioration. They can be delivered safely anywhere in the world, the end window types being sealed in aluminium containers, as illustrated on this page. They can be tested on reception (with gamma radiation) without breaking the seals.



While several more or less tubes in excess of the value of going to great lengths to ensure that the tubes and accessories going in the publication are accurate, the right is reserved to make any necessary changes to the details as we may consider to be in the interests of the customer.




**MULLARD LTD**  
X-RAY DIVISION  
New Road, Mitcham, Surrey  
England

Telephone: Mitcham 3411  
Telex: Mullard, London  
Mullard, New

Printed in England

MX 820



## Mullard halogen quenched Geiger- Müller tube Guide



# Reliability . . .

. . . largely depends on the careful selection and use of components. This folder is a guide to the most reliable Geiger-Müller tubes available.

Mullard halogen quenched tubes are in a class apart for strength, dependability and consistent quality. There is never any need to order matched sets, as tubes of any one type operate at the same working voltage.

They withstand reasonable over-voltage, incorrect polarity, and, in the case of all metal types, surprisingly severe vibration and drop tests.

They are very long lived, and may be stored for years without deterioration.

They can be delivered safely anywhere in the world, the end window types being sealed in aluminium containers, as illustrated on this page. They can be tested on reception (with gamma radiation) without breaking the seals.



Whilst every care is taken to ensure that at the time of going to press the illustrations and specifications given in this publication are accurate, the right is reserved to make improvements or to modify such details as we may consider to be to the advantage of the customer.



**MULLARD LTD**

X-RAY DIVISION

New Road . Mitcham Junction . Surrey  
England

Telephone:  
Mitcham 3471

Telegrams: Emvalco,  
Mitcham, Telex