

**Glycotrack : carbohydrate detection kit / Oxford GlycoSystems Ltd.**

**Contributors**

Oxford GlycoSystems Ltd.

**Publication/Creation**

Abingdon, Oxfordshire : Oxford GlycoSystems, [1995]

**Persistent URL**

<https://wellcomecollection.org/works/hjduhra6>

**License and attribution**


Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome  
collection**

Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

**GlycoTrack™**

*Carbohydrate Detection Kit*



**Cat. No:** K-050

**Price:** £95.00

- Application:** The positive detection of carbohydrate moieties in glycoproteins by specific labelling and membrane blotting.
- Sample Range:** From 10 ng to 20 µg of glycoprotein in 20 µl may be analyzed depending on the nature and extent of glycosylation in the glycoprotein.
- Number of Assays:** Reagents are sufficient for more than 50 labelling reactions and over 26 membranes each with an area of 1600 mm<sup>2</sup>.
- Starting Material:** Samples from cell lysates to pure proteins can be analyzed in conjunction with SDS-PAGE/ Western blotting or dot blotting.
- Storage:** The kit should be kept at 4°C before use. Components which must be kept at -20°C after reconstitution are grouped together in a separate foam insert for easy storage.
- Description:** A fundamental question asked by all scientists studying eukaryotic proteins is: "Is my protein glycosylated?" GlycoTrack has been specifically designed to provide an easy, definitive, answer to this question.
- Oxidation of the sample glycoprotein with periodate and subsequent reaction with biotin-hydrazide leads to the specific incorporation of biotin into the carbohydrate moiety. The biotinylated glycoprotein is then detected with a streptavidin-alkaline phosphatase conjugate. Two protocols for the detection of glycoproteins are provided:
- Protocol 1** has been designed for labelling glycoproteins already immobilized on a membrane and is particularly suitable when multiple samples (e.g. column fractions) are to be analyzed or interference by substances in the glycoprotein samples is suspected.
- Protocol 2** has been designed for solution labelling before electrophoresis and blotting and is the more sensitive method as all generated aldehydes are available for biotinylation.
- Kit Contents:**
- ❖ Periodate
  - ❖ Hydrazide
  - ❖ Blocking reagent
  - ❖ Protocol
  - ❖ Bisulphite
  - ❖ Streptavidin-alkaline phosphatase conjugate
  - ❖ 5-bromo, 4-chloro, 3-indolyl phosphate
  - ❖ Nitroblue tetrazolium
  - ❖ Native markers
  - ❖ Biotinylated markers
  - ❖ Ribonuclease B

- Specificity:** Minimal non-specific interactions through the specificity of periodate oxidation and biotinylation.
- High Sensitivity:** An extremely low background and high sensitivity is achieved through the Biotin-Streptavidin Detection System. Sample amounts from 10 ng to 20 µg of glycoprotein in 20 µl can be analyzed with as little as 1 ng of glycoprotein being detected.
- Convenient:** All essential reagents are included.
- Easy to Use:** Complete step by step protocols for both membrane and solution labelling are included together with a guide to interpreting results.
- Low Cost:** Reagent quantities and number of assays per kit have been optimized to minimize cost per analysis.
- Rapid:** Your question "Is my protein glycosylated?" can be answered in less than 5 hours.



**UK and Europe:** Oxford GlycoSystems Ltd, Unit 4, Hitching Court, Blacklands Way, Abingdon, Oxon OX14 1RG, UK.  
**Telephone:** +44 (0)235-553066 **Fax:** +44 (0)235-554701

**USA and Canada:** Oxford GlycoSystems, Inc., Cross Island Plaza, 133-33 Brookville Boulevard, Rosedale, NY 11422, USA.  
**Telephone:** (718) 712-2693 **Tollfree:** 1-800 722 2597 **Fax:** (718) 712-3364