Hetofrig CL410: -80C, 460 litres upright freezer is characterized by... / Heto-Holten A/S.

#### **Contributors**

Heto-Holten.

### **Publication/Creation**

Allerød: Heto-Holten, 1993.

#### **Persistent URL**

https://wellcomecollection.org/works/yrw7fhp3

#### License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection 183 Euston Road London NW1 2BE UK T +44 (0)20 7611 8722 E library@wellcomecollection.org https://wellcomecollection.org

## **HETOFRIG CL 410**

# -80 °C,460 litres Upright Freezer is characterized by:

Maximum storage volume: of 415 litres, total volume being 460 L on 0.75 m<sup>2</sup> of floor space only.

Convenient to use – easy access to products through a unique draw-shelf system.

Minimum maintenance – maximum energy saving thanks to a heavy insulation and an efficient refrigeration system.

Maximum safety: – simple to operate via the touch sensitive membrane panel. The freezer is equipped with the most advanced microprocessor control and alarm system

available today.
All commonly
required
parameters for
running
temperature and
alarm limits are
preset requiring
no further
adjustment
unless
specifically
needed.







## *HETOFRIG CL 410*

Control system: consists of two independent controllers a master and a back-up. The freezer temperature is constantly monitored and controlled, and a comprehensive self-check routine is carried out with short, regular intervals, checking the temperature sensors, compressors, mains power, battery voltage, CO2 vessel content and the function of the back-up controller.

When detecting a fault, the system immediately starts the alarm and, in case of a faulty controller, switches over to the operative controller. 12 possible sources of error are checked and diagnosed. A powerful self-charging

battery back-up ensures that all important information is available, even after 72 hours' power failure.


Alarm system: consists of an acoustic and a visual alarm. Connection to external alarm possible. Memory for several sets of alarms, each comprising type, date and time for the alarm, as well as the highest freezer temperature during each alarm state. The acoustic alarm can be interrupted by a password. This is particularly useful on freeze start up or

during long loading or unloading periods, when the door will be open, thus avoiding unnecessary alarm



An alpha-numeric data printer is available (cat. No. 868310). The printer enables full documentation of cabinet temperature, time, date, and control parameter settings. Time/date and diagnoses for alarms as well as max, temperature during

alarm are printed. A paper reeler for the data printer is also available (cat. No. 868311).



### **Touch Sensitive Membrane**


The freezer system is operated through a touch sensitive membrane

through a touch sensitive membrane panel with a smooth wipe surface. All temperature settings are protected by a password preventing unauthorized temperature changes. Set points for operating temperatures, alarm, CO<sub>2</sub> back-up and pre-freezing are built-in with

pre-set values, and give safety

#### Pre-freeze function:

The unique pre-freeze function raises the alarm and CO<sub>2</sub> back-up temperatures to a selected level for a limited, selectable period of time to avoid activating the alarm when loading large volumes of unfrozen



#### Test button:

A test button allows checking of the alarm systems and the CO<sub>2</sub> back-up without opening the freezer door as required by the FDA authorities.



#### Refrigeration system:

The refrigeration system is a cascade, air-cooled, hermetically sealed system. Powerful step-starting Danfoss compressors ensure fast start-up. The stepstarting function eliminates initial loading procedures, reduces power drain and extends compressor life. The system works in an excellent way even when the voltage variation

exceeds ±15%. It works at 30 °C and if ambient temperature constantly exceeds 32 °C, water-cooled refrigeration is recommended. Available as standard accessory



#### Door Sealing:


The door sealing is a unique seal with 3 flanges. The rubber seal remains soft – even at temperatures inferior to –90 °C, thus ensuring minimum loss of cold to the surroundings and avoiding ice crystal formation. Should the door seal happen to be damaged, it can easily be replaced.



The insulation of the freezer is 125 mm foam-in-place polyurethane on all sides, which means minimum uptake of heat from the surroundings, saving energy, extending the durability of the compressors, and offering maximum temperature stability (±0.5 °C). It gives minimum CO<sub>2</sub>

consumption during power failure and hence extends the time the products can be kept at a safe temperature. Normal CO<sub>2</sub> consumption at room temperature: Aprox 1-2 kg CO<sub>2</sub>/hour.

## 460 LITRES UPRIGHT FREEZER





Cryo storage system

A new generation of cryo storage racks is available for storage in the -80 °C freezers. The racks are easily assembled with strong extension nuts and can be disassembled when needed. This gives possibilities to:

have different cryo box sizes in the same cryo rack, giving optimum use

of the whole rack.

- change from e.g. 2" cryo box to 3" without the need of investing in other rack sizes.

The Heto cardboard or plastic boxes fit the racks perfectly, minimizing the dead volume in the racks and in the freezer. The most cost-effective solution available







#### One-hand operation:

The door lock is made of stainless steel and is operated with one hand. It ensures perfect door-sealing

and a key lock ensures safe storage of valuable products.



### CO<sub>2</sub> back-up:

CO2 back-up:
A CO2 back-up system is available
as an accessory (cat. No. 003360).
The CO2 back-up system is fully
controlled and monitored through
the control system. CO2 back-up
temperature level and flow
regulation are selected through the

Balance for CO<sub>2</sub> vessel is also available as an accessory (cat. No. 868320). It ensures precise control of CO<sub>2</sub> vessel content and max. safety for the products.





The shelves are constructed so that they can be used as fixed shelves. A simple lock system opens the door

The shelves can also be easily drawn more than halfway out of the freezer thanks to a rail system enabling easy access even to the products stored on the rear of the

This makes loading and unloading

of sample very easy and saves a lot of time and trouble in routine use. The shelves are made of stainless steel (AISI 316). The four shelves are standard equipment.



The automatic vacuum release system quickly equalizes the pressure in the freezer with the surroundings giving an easy opening of the door



#### 80°C Freezers

Heto offer a wide range of -80°C Freezers in upright and chest-freezer versions. The chest-freezer range comprises volumes from 123 to 420 litres. All freezers are equipped with the advanced control and alarm system. —40 °C chest freezers are also available. (Please contact us for further information).

# **HETOFRIG CL 410**

Technical specifications	Model:	CL 410	
Cat. No.		003405	003355
Temperature range,	°C	-85 to -45	
Temperature stability,	°C	±0.5	
Uniformity,	°C	±1	
Storage volume/Total volume,	L	415/460	
External dimensions D×W×H,	mm	935×800×1970	
Internal dimensions D×W×H,	mm	675×540×1500	
Shelf dimensions D×W×H, mm		615×510×300	
Refrigeration		2 Danfoss Compressors	s SC 18 C, Cascade
Refrigerant		R502/R170	
Internal material		stainless steel (AISI 316)	
External material		painted steel	
Insulation		125 mm polyurethane foam on all sides	
Power consumption, average,	W	600	600
Power consumption, max., W		1300	1300
Power,	V/Hz	115/60	220/50
Weight,	Kg	294	

NB! The freezer is as standard delivered with a right-hand door. But the construction of the freezer enables an easy change from right-hand to left-hand door, which gives you a greater flexibility when deciding where to install the freezer.

#### Accessories

Cat. No.	Description	
868310	Alpha-numeric data printer	
868311	Paper reeler for printer	
868312	Paper roll for printer, one pc.	
868313	Ink ribbon for printer	
003360	CO <sub>2</sub> connection kit incl. all connections and 12 V battery back-up	
868320	Balance for CO <sub>2</sub> vessel	
003361	CL-flex-0 for upright freezer CL 410. Includes hand grip, bottom and extension clips Dim: D137×W137×H285 mm	
003362	One shelf for cryo storage racks, excl. clips Dim: D137 × W137 mm	
336106	Extension clips for 10 shelves, made of stainles steel	
801040	CLB-52*/5 cm cardboard cryo box without insert Dim: D135×W135×H50 mm	
801041	CLB-73*/7 cm cardboard cryo box without insert DIM: D135×W135×H70 mm	

Cat. No.	Description	
801042	CLB-10 4"/10 cm cardboard cryo box without insert Dim: D135 × W135 × H100 mm	
822565	CLB 149 cardboard insert for 49 cryo tubes	
822544	CLB I 64 cardboard insert for 64 cryo tubes	
822543	CLB I 81 cardboard insert for 81 cryo tubes	
822542	CLB I 100 cardboard insert for 100 cryo tubes	
822541	CLP-5 2"/5 cm plastic cryo box for 100 cryo tubes Dim: D135 × W135 × H50 mm	
822555	CLP-73"/7 cm plastic cryo box for 100 cryo tubes Dim: D135xW135xH70 mm	
822535	PRS 13 test tube rack for 13 mm tubes, D102 × W102 mr	
822536	PRS 16 test tube rack for 16 mm tubes, D126 × W126 m	
822537	PRS 20 test tube rack for 20 mm tubes, D128 × W103 mm	
822539	PRS 25 test tube rack for 25 mm tubes, D122 × W122 mm	
822556	PRS 30 test tube rack for 30 mm tubes, D109 × W109 mm	



Heto Lab Equipment A/S Gydevang 17-19 DK-3450 Allered Denmark Telephone: +45 48 14 27 77 Telefax: +45 42 27 46 55 Telex: 40 037