'Wellcome' brucella milk ring test antigen.

Contributors

Burroughs Wellcome and Company

Publication/Creation

1967.

Persistent URL

<https://wellcomecollection.org/works/cxv6r6f7>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

'WELLCOME' BRAND

BRUCELLA MILK RING TEST ANTIGEN

'Wellcome' Brucella Milk Ring Test Antigen is designed for use in the detection of *Brucella abortus* antibodies in milk samples. It consists of a suspension of killed *Br. abortus*, Strain 99, stained with hæmatoxylin.

METHOD OF USE

The milk samples should be tested in narrow glass tubes which will contain 1 ml. in a column about 2 cm. high. Standard tubes 2 in. or 3 in. long and $\frac{3}{8}$ in. diameter are suitable.

One drop (0.04 ml.) of 'Wellcome' Brucella Milk Ring Test Antigen is added to 1 ml. milk and gently mixed by inverting the tube several times. Mixing should be carried out within one minute of adding the Antigen. The presence of darker stained areas in the top layer of milk after standing for a few minutes indicates inadequate mixing.

The tubes are incubated for one hour at 37° C. and the results are interpreted as follows.

INTERPRETATION

Positive reaction : A definite blue cream ring with a white or slightly blue milk column.

Negative reaction : A cream layer lighter than or the same colour as the milk column.

ACCURACY

The milk ring test is mainly of use for screening herds by testing samples of bulked milk. Unreliable results are obtained from individual cows at the beginning or end of lactation and from cows affected with mastitis. Cows which have been vaccinated as adults with *Br. abortus*, Strain 19, usually give positive results.

STORAGE

The Antigen should be stored at 2° to 10° C. It should not be frozen.

PRESENTATION

'Wellcome' BRAND Brucella Milk Ring Test Antigen is supplied in bottles of 4 ml. (100 tests) with droppers.

Also available :

'Wellcome' BRAND Brucella Abortus Agglutinating Serum. A drop added to 1 ml. milk may be used to give a positive control in the ring test.

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(The Wellcome Foundation Ltd.)