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LEIPZIG.

DESCRIPTION

The colour plates are used for detecting defects of colour-vision. Disturbances of Red-Green-Perception are common: red and green signs are very important as they are used for Railway, Naval, Aircraft and Traffic Signals. Accordingly the greater number of the plates (Nos. 3-31) are intended for detecting defects of Red-Green-Perception. Tables 32-35 are for detecting the less common defects of Yellow-Blue-Perception.

A. Testing for Red-Green-Perception.

Plates 3-31

As it is important to detect defects of Red-Green-Perception, and also minor degrees of these, it is pointed out that tables 3-31 are designed to achieve this aim in different but supplementary ways. Therefore all the plates should be shown to the subject, so that assessment of his colour vision can be complete.

Plates 3-18

These show characters (Numbers or Letters) which are composed of variegated dots and a background consisting of different coloured dots. For the normal person the colours of the prominent detail contrast sharply with the background and they can be read by him without difficulty. (Normal Trichromatic Vision).

Those persons, however, whose red-green perception is defective will either perceive imperfectly, or not at all, the red and green tones which preponderate in the colour-blending of the dots; but the yellow, bluish and grey tones are quite visible to them. As these are the colours of the dots of which the background is composed, the difference between the detail and background is indistinct or absent for Red-Green-Defectives and they have difficulty in recognising the detail. Complete failure to recognise the detail is a result of Red-Green-Blindness; incorrect interpretation of the detail or the recognition of it after prolonged effort, intense gazing or oblique peering, shows Red-Green weakness. (Anomalous Trichromatic Vision).

Not all those suffering from Anomalous Trichromatic Vision have the same type of disturbance of colour-perception, therefore, some of the subjects will be able to recognise details which others cannot. For this reason it was necessary to assemble numerous types of plates, selected from the examination of a large number of colour-blind persons and those with weak-colour-vision.

It is pointed out that tables 13-18 show whether the subject is suffering from Protanopy (Protanomalie), i.e. red-green defect present, but principally a defect of perception of red; or from Deuteranopy (Deuteranomalie) i.e. red-green defect present but principally a defect of perception of green.

Plates are read as follows:-

<u>Protanopes</u>	<u>Deuteranopes</u>	<u>Normal Persons</u>
Plate 13 mostly not	nearly always	as 8
Plate 14 often	nearly never	as 65
Plate 15 mostly not	mostly	as 5
Plate 16 mostly	mostly not	as 6
Plate 17 as 4	as 3	as 34
Plate 18 as R	as Z	as RZ

Plates 19-22

These are based on the strongly marked powers of perception of differences in intensity of colour found in colour-blind subjects. With these plates, normal subjects distinguish without difficulty the details by their colouring from the background. These colour variations escape those with defective colour vision but, owing to their great powers of perception of differences of intensity they appreciate some lighter detail which appears hidden to ordinary subjects. In these tests a lot depends on the statement made at first sight by the subject.

These will be easily visible:-

<u>For normal subjects</u>	<u>For defective subjects</u>
Plate 19 letters OH	31
Plate 20 3	8
Plate 21 8	2
Plate 22 mostly no signs	16

Subjects who are severely colour-blind will recognise the very dark tones on a light-coloured background, not by the colouring but by the difference in intensity. On plates 1 and 2 the detail is very dark in contrast with the background and so they will be read correctly by everyone.

including the severely colour-blind. As these plates, on account of their easy perception by all subjects, are suitable for introduction in the use of the tests they are placed at the beginning. (For their use as simulation tests turn to page 9)

Plates 23-28

These are a supplementary test; they contain no numbers or letters, but only dots, in colours that are easily confused, in different shades and degrees of intensity. One dot is chosen and shown to the subject who has to pick out all equally coloured dots and to indicate them on the plate with a pointer. Normal subjects can do this very quickly but others will make mistakes, especially if it is ensured that the smallest dots of the chosen colour are also indicated. These plates serve chiefly for tracing defects of Red-Green-Perception. They have not only a qualitative but also, to some extent, a quantitative character.

Plates 29-31

These are based on a completely different principle from those hitherto discussed. They make use of the contrast-effect which is increased in defective people, to detect Red-Green defects. Many of them state that the grey dots, surrounded by a light-green surface, as in Plate 29, give them the same impression as if they were reddish, reddish-grey or violet dots. They, therefore identify the grey dots on Plate 29 as reddish, reddish-grey or violet.

This contrast effect can be seen more often and more clearly if one shows brown or brownish dots on a bright red background (Plate 30). To defective people they seem to be green, grey-green or olive coloured. The figure 5 on Plate 31 is often identified as greenish. This test can be very important for defects of Red-Green-Perception especially if, while using other Plates, suspicion has already been aroused. During the test the subject should be asked to state the shade of colour of the central dot of the green (Plate 29) or of the red (Plate 30) squares. The figure 5 of Table 31 can easily be recognised by the subject but it is important to know in which colour he sees it.

TESTING OF BLUE-YELLOW-PERCEPTION.

Plates 32-35

These set out to detect disturbances of Blue-Yellow-Perception. This is very uncommon but among those suffering from Red-Green Disturbance there are sometimes subjects whose perception for Blue-Yellow is also defective. In other rarer cases Red-Green-Perception can be quite normal and Blue-Yellow-Perception defective. The tracing of such cases is not only of theoretical interest. Those who suffer from such defects cannot distinguish blue-green or blue; yellow or orange. This might become obvious in occupations with exacting demands, e.g. railwaymen. It is important, therefore, that our colour-vision tests should include also Plates for examining blue-yellow

perception, (in contrast with other types of pigment tests). These plates should be shown after the completion of the Red-Green-Perception test, because without knowing the condition of the latter, a judgement on the blue-yellow sense would be incomplete.

Plates 32 and 34 will be read by normal subjects without any effort, but subjects with more serious defects of Blue-Yellow-Perception cannot read them at all or only imperfectly. Tables 32 and 35 are intended to detect less seriously defective Blue-Yellow-Perception.

Plate 33 shows on an orange background the numbers 4 and 5 in pink, which can be recognised by a normal person at about 1 metre. Among these there are in lighter dots the Numbers 9 and 2, which are, however, not visible to a normal subject. Those whose blue-yellow vision is defective (~~Tritanopia~~) read quickly these numbers but cannot decipher the figures 4 and 5. Those who are able to distinguish colours can see on Plate 35 the figures 4 and 9 in blue dots. Among these is shown a 2 in lighter colouring which can be recognised by normal subjects with some effort. Tritanomalous subjects read only the 2; but cannot see the figures 4 and 9.

II DIRECTIONS FOR USE.

Generally the Plates should be shown at a distance of 1 metre; on account of the smallness of the dots on Plates 23-28 the distance can be shortened to 60 cm. The best lighting is diffuse daylight, not sunlight.

Artificial light, especially the so-called electric daylight lamps, are also suitable for Plates 1-18. Plates 19-35 give different results when they are shown in artificial light instead of daylight. Like the normal subjects the examinee should look at the Plates and then say quickly what he sees.

Prolonged effort, moving the head or drawing the digits with the finger in the air, signify defects of colour-vision.

The order in which the plates are shown can be chosen as required, but it is recommended that plates 1 and 2 be shown at the beginning, as they can easily be read by colour-blind persons and are therefore suitable as an introduction to the test.

Plates, showing figures or letters, permit an extensive examination. The examiner is, however, not limited by mere identification of digits, but can also get the subject to point out similarly coloured dots, especially Plates 23-38. If doubts arise the subject can be asked to draw the figures or letters with a pointer. A normal person can do this correctly but others become confused by the colour mixture of the background and lose the digit. Plates 2, 11 and 12 are suitable for drawing the digit. It is wisest to show the Plates in a different order to each subject, so that they cannot be guessed.

The mistakes observed should be noted, but should not be told to the subject. If there are many errors, i.e. digits are not recognised, it will be obvious that there is defect of colour vision. If the errors are

unimportant, e.g. similar numbers like 3 and 5, 6 and 9 are confused, these Plates should be noted and shown again at the end of the examination. If these errors are repeated these slight variations indicate some defect.

III DISCLOSING OF WRONG STATEMENTS OF THE EXAMINEE

A. Simulation

If the examinee states that he cannot recognise any of the tests and that he is quite colour blind, he will, in most cases, also be suffering from some other eye disease; such as optical weakness, Nystagmus, Photophobia.

In the absence of these conditions, there are grounds for suspecting simulation of colour-blindness. Simulation will be proved, if the subject does not interpret Plates 1 and 2 correctly; these are visible to all on account of the differences in intensity.

Malingers mostly state that they cannot recognise red or green. To detect simulation, any Plates are shown to the examinee and his statements are noted. Tables 13-18 chiefly confuse malingerers as they are read differently by Protanopes and Deuteranopes (described on page 3).

There are reasonable grounds for suspecting simulation if the subject maintains that he cannot read Table 34 with the yellow figures on a blue background, as most people suffering from red-green defects read this plate without difficulty. Proof of simulation is failure to read Plates 1 and 2 which can be read even by red-green-blind subjects owing to the differences

in intensity. Tables 19-22 can be used also as real red-green defectives can read the figures.


B. DISSIMULATION.



It often happens that a subject tries to conceal his colour-blindness. The memorising which might be possible with a small selection of Plates, is not possible owing to the great number and variety of these Plates. To avoid a fortunate guess, Plates should not be shown in order. They can also be turned upside down and with many plates quite readable details can be seen in this way, but memorising the plates is useless.



Sketching the letters and figures, as mentioned above is a good test of the state of colour vision, especially Plates 2, 11 and 12.



Finally one can change the Plates by colour filters (e.g. glass plates). Some Plates (Numbers 9, 10, 12, 14 and 18) are quite unreadable by normal persons through a blue or red colour filter. If the figures are identified under these conditions the subject has learnt them by heart. Other filters in pink or green cause several variations so that the examiner may catch out such a subject more easily (Filter trap).














23

24

25

26

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This work has been identified as being free of known restrictions under copyright law, including all related and neighbouring rights and is being made available under the Creative Commons, Public Domain Mark.
Copyright 1939 by GEORG THIEME, LEIPZIG