

**Trojan high frequency (violet-ray) electrical outfits for health and beauty /
John Shaw & Sons, Wolverhampton, Ltd.**

Contributors

John Shaw & Sons (Wolverhampton, England)

Publication/Creation

Wolverhampton : John Shaw & Sons, [between 1910 and 1919?]

Persistent URL

<https://wellcomecollection.org/works/ycapf7ky>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome
collection**

Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>



HIGH FREQUENCY
(VIOLET-RAY)

Electrical Outfits

— for —

HEALTH AND BEAUTY.

INTRODUCTION.

HIGH Frequency Electric Current, very generally known as Violet-Rays, was discovered towards the end of last century, and to-day its remarkable curative powers are universally acknowledged. It must not, however, be confused with the Ultra-Violet or Sun-Ray treatment, which is entirely different and ought only to be indulged in under the direction of a qualified specialist.

High Frequency is the collective expression denoting alternating electric currents, the direction of which change with great rapidity—equal to more than a million alternations per second. This allows the actual application of electricity—with all its curative and sedative powers—to the nerves of the body without shock or any other irritating effect, bringing almost instant relief to the innumerable ailments to which the human frame is subject.

These currents or Violet Rays are applied by means of a special apparatus—the “TROJAN” Outfits herein described having been designed—not only to provide the medical profession with a complete and dependable outfit, but also with such simplicity that they can be easily and simply operated by an entirely inexperienced person to treat his or her own complaint personally. Full and detailed instructions derived after careful medical study for treating a wide range of ailments and for generally promoting health and beauty follow hereafter.

Highly recommended and can be indulged in with absolute safety and confidence.

John Shaw & Sons, Wolverhampton, Ltd.

ESTABLISHED 1795.

Head Office and Works:
WOLVERHAMPTON.

Phone: 1001.

Grams: “Shaw,” Wolverhampton.

London Office and Showrooms:
29/30, SHOE LANE, E.C. 4.

Phone No.: Central 2163.


Discount 25%



High Frequency Electricity

— for —

BEAUTY and HEALTH.



BEAUTY.


BEAUTY and Hair Specialists strongly recommend the "TROJAN" High Frequency Outfits. This treatment creates health and vitality in the skin, removes wrinkles, creases and blemishes from the face, stimulates the circulation of the blood and makes the flesh firm, creating that healthy glow which is the foundation of beauty.

High Frequency treatment for the hair has wonderful effects; dandruff and falling hair are quickly cured and the growth stimulated.

HEALTH.

THE "TROJAN" High Frequency treatment is strongly recommended for all nervous troubles; pains such as headaches, neuralgia, chilblains, etc., are instantly relieved after the first treatment—such chronic diseases as rheumatism, lumbago, sciatica, catarrh, nephritis, etc., etc., are completely cured after a course of treatment.

High Frequency Electricity is Life!




INSTRUCTIONS.

FIRST plug into the nearest electric lamp holder. Select the electrode recommended for the complaint and push right home into the operating handle and switch on the current. Next turn the dial knob until a violet radiance appears in the Electrode. The machine is now ready for operation. A further turn of the dial knob will increase the strength of the current as desired. When treatment is finished turn the dial knob backwards to zero before switching off. Never fit or remove the Electrode without first turning off the current, i.e., without first turning the dial knob back to zero.

TROJAN

High Frequency Electricity




MODEL No. A.—Complete in strong wood carrying case, covered with black leatherette, plush lined, fitted with four superfine velvet finish glass electrodes and metal saturator, operates on any voltage from 110—250.

Price - £2 10s. 0d. complete.

MODEL No. D.—Complete in strong wood carrying case, covered with crocodile leatherette, plush lined, fitted with nine superfine velvet finish glass electrodes, metal saturator, ozone generator and test tube. Operates on any voltage from 110—250.

Price - £4 17s. 6d. complete.




The nickel plated clip shown in this set is for attaching to wrist to prevent interference with wireless reception.




High Frequency Electrodes.


THE whole of the electrodes illustrated below are included in the "TROJAN" Model "D" High Frequency Outfit, and in the Model "A" the following numbers 0, 3, 4, 10, 104, only are supplied. A treatment chart and complete instructions for using the various Electrodes is included in every outfit.




No. 0.—Surface Electrode, for general application to any part of the surface of the body or face.




No. 2.—Intensifying Electrode, similar to No. 0, but for treating the more deep seated complaints, such as Rheumatism, etc.




No. 4.—Comb Electrode, for promoting the growth of the hair and general scalp treatment.




No. 3.—Throat Electrode, specially shaped for external application to the neck and throat.




No. 10.—Fulguration Electrode. This Electrode is for treating warts, moles, etc.




No. 30.—Spiral Electrode. Specially designed for the treatment of heart diseases.



No. 35.—Small Concave Surface Electrode, for treating protruding parts, such as the chin, knuckles, etc., and can also be used for applying chemical solutions in conjunction with High Frequency.




No. 27.—Rectal Dilator Electrode can be supplied in three diameters, viz.: 2-in., 1-in., 1/2-in.




No. 104.—Metal Saturator Electrode, to be held by patient during massage.



No. 102.—Glass Roller Electrode, for High Frequency massage.



No. 105.—Test Tube. The object of this Tube is to verify that the High Frequency current is passing through the whole of the patient's body during massage. The patient holds the No. 104 Electrode and the operator passes the Test Tube along any part of the body (not actually in contact) and it will be noted that this Test Tube will become luminous, thus indicating that the current is passing from the Electrode right through the entire body.



No. 101.—Ozone Generator Electrode. As its name implies, this Electrode is for generating ozone in conjunction with High Frequency. A few drops of "TROJAN" Ozone Inhalant, as described below, are put on the cotton wool in the glass bulb and the mixture generated is inhaled in the form of vapour directly into the lungs through the nostrils, destroying all germs and generally soothing any inflamed membranes. Supplied complete with rubber tube and bulb—which latter when regularly pressed causes a continuous flow of Ozone Vapour to the nostrils. The operator must hold the glass stem of Generator when under treatment.



"TROJAN" High Frequency Lubricant, for softening the surface of the skin and generally lubricating, allowing the Electrodes to move over the part treated with ease.
Price .. 2s. 6d. per bottle.



"TROJAN" Ozone Inhalant, for use in conjunction with the Ozone Generator described above. Price 1s. 9d. per bottle.