

**"Burdizzo" pincers for the bloodless castration of sheep, bulls, pigs, horses, dogs, etc. : recommended by the veterinary surgeons ... : inventor and sole maker: Dr. N. Burdizzo, Turin (Italy).**

**Contributors**

N. Burdizzo (Firm)  
Arnold & Sons.

**Publication/Creation**

Turin? : N. Burdizzo, [ca.1930]

**Persistent URL**

<https://wellcomecollection.org/works/rydsf5pd>

**License and attribution**

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).




Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

TRADE **“BURDIZZO,”** MARK

PINCERS  
 FOR THE BLOODLESS CASTRATION OF  
**SHEEP, BULLS, PIGS, HORSES, DOGS, Etc.**  
 RECOMMENDED BY THE VETERINARY SURGEONS

NICKEL  
 PLATED.

FOUR  
 DOUBLE  
 LEVERS.



VERY  
 SOLID  
 CONSTRUCTION.

CAN BE CARRIED  
 IN THE POCKET.

NOT GENUINE and  
 Not RELIABLE without  
**Dr. BURDIZZO'S**  
 (INVENTOR AND SOLE MAKER)  
 TRADE MARKS:



N. BURDIZZO  
 LA MORRA (ITALIA)

GUARANTEED TO BE OF PERFECT MANUFACTURE  
 AND TO GIVE SATISFACTION.

- |  |                |
|--|----------------|
| <b>BABY «BURDIZZO»</b> (about 8in. long, 1lb. weight), for Lambs, from 1 week to six months old  | <b>£ 2.5.0</b> |
| <b>SMALL SIZE</b> (about 12in. long 3 lb. 4 oz. weight), for lambs, for rams from 3 months to any age, for Bulls under 9 months, for Pigs, Dogs                      | <b>„ 3.0.0</b> |
| <b>LARGE SIZE</b> (about 13 inches long, 4 lbs. 6 oz. weight) for all classes of grown animals, but particularly for Horses and Bulls of any age                     | <b>„ 3.0.0</b> |
| <b>ELONGATED SIZE</b> (about 16in. long, 4 lbs. 8 oz. weight). This is the «large size» with elongated handles. Particularly useful in bulls with thick cord and sac | <b>„ 3.0.0</b> |

**POSTAGE EXTRA.**

*Pincers may be sent C.O.D. by Post.*

**EVERY FARMER CAN CASTRATE HIS BULLS, RAMS,  
 ETC. WITH ONE ASSISTANT, AS BULL CALVES  
 ARE OPERATED STANDING.**

- Certain** — The result is invariably a success.
- Convenient** — The castrated animals can be sent directly to pasture without fear of after-effects or need of further treatment.
- Rapid** — Bulls and Rams can be castrated in one or two minutes. With sufficient assistants 50 Bull calves can be castrated per hour with one instrument.
- Humane** — There is no bleeding, the bag is not cut, the cords are crushed under the skin.
- Safe** — Animals can be castrated with safety whether the day be rainy or windy and dusty, or the pasture be full of flies.
- Economical** — The same instrument will castrate every kind of domestic animals whatsoever its age.
- It avoids screw-worms worries** to farmers. Castration can be performed all year round in Blow-fly infested farms.
- It brings home one pound sterling profit per head** in slaughter stock.
- Young Lambs are easily castrated and tailed without bleeding.**

**INVENTOR AND SOLE MAKER: Dr. N. BURDIZZO - TURIN (Italy).**

DISTRIBUTORS

**ARNOLD & SONS**  
 50-52, WIGMORE STREET,  
 LONDON, W. 1.

Printed in Italy.

# CASTRATION and DOCKING LAMBS in NEW ZEALAND

with The **BURDIZZO** "Wonderful Success,,  
THE PICTURE: THOMSON BROS., KAIRANGA,  
NEW ZEALAND, IN OPERATION.

NO LOSS  
OF BLOOD

—  
NO WOUND

—  
CLEAN  
AND SORE



THE  
TESTICLES  
WILL  
WITHER  
AWAY  
IN THE SACK

THE WAY: CLOSE PINNERS ON SACK.  
OPEN AND LET THE LAMB RUN.

## "2500 Lambs without the loss of a single one"

KAIRANGA, N.Z., Sept. 15, 1927.

Messrs. Abraham & Williams Ltd.,

Agents «Burdizzo», Palmerston North.

*We have to advise having used the Burdizzo Pincers (small size) for two seasons, with wonderful success. Last season we docked 2,500 ram lambs without the loss of a single one.*

*We consider that this new method of docking is far superior to the old method. It is completely sanitary, and avoids all chance of blood poisoning. The lambs do not get such a set-back, and will feed almost right away after the operation.*

*One feature of the new method is the possibility of docking in any weather, wet or fine, which is not advisable with the knife.*

*This method of docking is much faster than with the knife.*

(Sgd.) Messrs. THOMSON BROS.

## A Few Statements from our 16 pages Booklet:

F. D. McMASTER, DALKEITH, CASSILIS, N. S. W.

*«I have personally given them an exhaustive trial and find them satisfactory in every way and all you claim for them.»*

WALTER HAWKER, ANAMA STATION, CLARE, S. A.

*«Please send by parcels post one pair of Burdizzo Pincers. I have used a pair at Kalabity with very good results.»*

From A. D. MACKINNON, DALNESS, EVANDALE, Tas.

*«I have given them a good trial on rams of all ages and find them quite satisfactory. My losses have been none though I used them in the middle of summer.»*

From R. FALCONER, YULGAN, BOLGART, W. A.

*«Have operated on 900 ram lambs this year with entirely satisfactory results. We would not go back to the knife on any account.»*

## Docking a lamb's tail by the Burdizzo method:

*Twenty-eight hundred lambs were docked by the Burdizzo method by the Spring Valley Water Company, San Mateo County, California, U.S.A., during month of February, 1929, with splendid success.*