

(H2O)[to the power of infinity] = RX : True or false? : TRUE! Thanks to ELIX continuous deionisation / Millipore.

Contributors

Millipore Corporation.

Publication/Creation

[Bedford, Mass.] : Millipore, 1994.

Persistent URL

<https://wellcomecollection.org/works/eczugyft>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

Put our claim to the test !

Return the enclosed reply-paid card to receive full information on the Milli-RX Analytical Grade Water Purification System - and don't forget to ask for our new poster!



If you would like to receive our new 42cm x 60cm poster, please tick the appropriate box on the reply-paid card.

To obtain additional information, please use the enclosed reply-paid card or call your nearest Millipore office :

AUSTRIA, CENTRAL EUROPE, C.I.S., AFRICA, MIDDLE EAST and GULF
Millipore Ges.m.b.H.
Hietzinger Hauptstraße 145
A-1130 Wien, Austria
Tel. (43) 1 877 8926
Telex 131464 miliv a
FAX (1) 877 1654

BALTIC REPUBLICS
Millipore Oy
Lantkatu 2
FIN-02770 Espoo, Finland
Tel. (90) 8045 110
FAX (90) 8596 616

BELGIUM and LUXEMBOURG
Millipore S.A.-N.V.
Rue de la Fusée 60/
Raketstraat 60
B-1130 Bruxelles-Brussel
Tel. (02) 242.17.40
FAX (02) 242.27.85

CZECH REPUBLIC
Millipore (Austria) Ges.m.b.H.
Obchodní Zastoupení
Ricanova 21
CZ-169 00 Praha 6
Tel. (42) 2-35 02 27
(42) 2-35 23 75
FAX (42) 2-35 02 27
(42) 2-35 23 75

DENMARK
Millipore A/S
Baldersbuen 46
DK-2640 Hedehusene
Tel. 46 59 00 23
FAX 46 59 13 14

FINLAND
Millipore Oy
Lantkatu 2
FIN-02770 Espoo
Tel. (90) 8045 110
FAX (90) 8596 616

FRANCE
Millipore S.A.
BP 307 - F-78054 Saint-Quentin
Yvelines Cedex
Tel. (1) 30.12.70.00
Telex 698 371 F
FAX (1) 30.12.71.80

GERMANY
Millipore GmbH
Postfach 5647
D-65731 Eschborn
Tel. (06196) 494-0
Telex 40 72 704 mil d
FAX (06196) 482237

HUNGARY
Millipore kft
Andor u 47-49
H-1119 Budapest
Tel. (1) 166 86 74
Telex 22 30 57
FAX (1) 185 32 35

ITALY
Millipore S.p.A.
Via B. Buozzi 4
I-20090 Vimodrone MI
Tel. (02) 25078.1
FAX (02) 2650324

THE NETHERLANDS
Millipore B.V.
Penningweg 33
Postbus 166
NL-4870 AD Etten-Leur (N.B.)
Tel. (01608) 22000
Telex 54 840 milwa nl
FAX (01608) 22436

NORWAY
Millipore AS
Postboks 39 Simensbråten
N-1107 Oslo
Tel. 22 67 82 53
FAX 22 68 53 15

POLAND
Millipore (Austria) Ges.m.b.H.
Oddział w Warszawie
Ul. Jasnowarska 7
PL-01-745 Warszawa
Tel. (48) 2 669 12 25
(48) 2 663 70 31
(48) 2 663 70 32
Telex 825 422
FAX (48) 2 663 70 33

SPAIN and PORTUGAL
Millipore Ibérica, S.A.
Avda. del Llano Castellano, 13-3º
E-28034 Madrid
Tel. 91-729 03 00
FAX 91-729 29 09

SWEDEN
Millipore AB
Box 7090
S-172 07 Sundbyberg
Tel. 08-628 69 60
FAX 08-628 64 57

SWITZERLAND
Millipore AG
Chriesbaumstrasse 6
CH-8604 Volketswil
Tel. (01) 945-3242
Telex 82 94 08 MILAG CH
FAX (01) 945-3766

U.K. and IRELAND
Millipore (U.K.) Ltd.
"The Boulevard"
Blackmoor Lane
Watford, Hertfordshire WD1 8YW
Tel. (0923) 816375
Telex 24191 milpor g
FAX (0923) 818297


Cat. No. DE16/U. Printed in France, 08/94. Copyright © 1994, Millipore Corporation, USA.
Millipore is a registered trademark of Millipore Corporation. ELIX and Milli-RX are trademarks of Millipore Corporation.

$(H_2O)^\infty = RX$

True or False?

TRUE!

Thanks to ELIX™ Continuous Deionisation



If you need a more economical, environmentally friendly system which produces Analytical Grade Water from ordinary tap water, there's only one answer:

The new Milli-RX™ System from Millipore.

Unlike conventional water purifiers, the Milli-RX System never requires time-consuming regeneration using aggressive chemicals. Neither does it call for constant maintenance like distillation-based equipment.

Instead, Milli-RX combines proven Reverse Osmosis with our unique ELIX technology which uses an electric current to regenerate the system continuously. So it's clean, efficient and safe!

Every day, the Milli-RX system produces up to 1000 litres of pure water* from potable tap water:

- Resistivity: ≤ 15 Megohm·cm
- T.O.C.: ≤ 50 ppb
- Bacteria: < 5 cfu/ml

* Analytical Grade Water - Grade 2 ISO 3696/BS3978

This water quality is suitable for a wide range of applications:

Typical Applications

- ✓ Purified Water (European Pharmacopoeia)
- ✓ Microbiological Media Preparation
- ✓ Buffer and Analytical Reagent Preparation
- ✓ Pure Water Supply to Instrumentation

The Milli-RX system is simple to use, and the quality of water produced is constantly monitored.

No other water purification system can equal the new Milli-RX for performance, cleanliness and economy.

Please return the enclosed reply-paid card to obtain full information on the new Milli-RX water purification system.

MILLIPORE