The solution to cost effective sample preparation just became more attractive! / Waters Corporation.

#### Contributors

Waters Corporation.

#### **Publication/Creation**

[Milford, Mass.] : Waters Corporation, 1994.

#### **Persistent URL**

https://wellcomecollection.org/works/eaaa3cag

#### License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection 183 Euston Road London NW1 2BE UK T +44 (0)20 7611 8722 E library@wellcomecollection.org https://wellcomecollection.org

# The solution to cost effective sample preparation just became more attractive !

#### Austria and European Export (Africa, Central Europe, CIS and Middle East)

Waters Ges.m.b.h. Hietzinger Hauptstrasse 145 A-1130 Wien Tel. 1 - 8771807 Fax 1 - 8771808

#### Belgium and Luxemburg

Waters S.A. - N.V. Raketstraat 60 Rue de la Fusée 60 1130 Brussels Tel. 02 - 7261000 Fax 02 - 7261100

#### CIS

Waters Representation Office Ul. Miklukho-Miklaia 16/10 117871 Moscow Tel. 7 - 095 - 3307447 Fax 7 - 095 - 3107007

#### Czech Republic and Slovakia

Waters Ges.m.b.H. Ricanova 19 16900 Praha 6 Tel. 42 - 2 - 35 23 75 Fax 42 - 2 - 35 02 27

#### Denmark

Waters A/S Baldersbuen 46 2640 Hedehusene Tel. 46 - 598080 Fax 46 - 598585

#### Finland

Waters Oy Kuparitie 1 00440 Helsinki Tel. 358 - 0 - 5064140 Fax 358 - 0 - 50641411

#### France

Waters S.A. B.P. 608 78056 St. Quentin en Yvelines Cedex Tel. 1 - 30487200 Fax 1 - 30487201

#### Germany

Waters GmbH Hauptstraße 87 D-65760 Eschborn Tel. 06196 - 400600 Fax 06196 - 482388

#### Hungary

Waters kft Andor ut 47-49 1119 Budapest Tel. 36 - 11620686 Fax 36 - 11853236

#### India

Waters Pvt. Ltd. 1st. Floor Shankar House No. 1 R.M.V. Extension Mekhri Circle Bangalore 560080 Tel. 91 - 80 - 3341944 Fax 91 - 80 - 3340527

#### Italy

Waters S.r.l. Via B. Buozzi 4 20090 Vimodrone MI Tel. (02) 250781 Fax (02) 2650324

#### Norway

Waters AS Enebakkveien 119 0680 Oslo 6 Tel. 22 - 195360 Fax 22 - 195362

#### Poland

Waters Sp.zo.o ul. Jasnodworska 7 01 745 Warszawa Tel. 48 - 26691225 Fax 48 - 26637033

#### Spair

Waters Cromatografia S.A. Entenza 24 - planta baja 08015 Barcelona Tel. 93 - 3259616 Fax 93 - 3259896

#### Sweden

Waters Sverige AB Box 7024 172 07 Sundbyberg Tel. 08 - 6282440 Fax 08 - 6282455

#### Switzerland

Waters AG Chriesbaumstrasse 6 8604 Volketswil Tel. 01 - 9462666 Fax 01 - 9453790

#### The Netherlands

Waters Chromatography B.V. Postbus 379 Penningweg 33 4870 AJ Etten-Leur Tel. 01608 - 22000 Fax 01608 - 22436

#### UK and Ireland

Waters (UK) Ltd. The Boulevard Blackmoor Lane Watford, Herts, WD1 8YW Tel. 01923 - 816700 Fax 01923 - 219012

#### All other countries

Waters Corporation 34 Maple Street Milford MA 01757 USA Tel. (508) 478 - 2000 Fax (508) 478 - 6963

Sep-Pak and Waters are trademarks of Waters Corporation, USA.

Printed in France, 11/94. © Copyright 1994 Waters Corporation, USA. P949405





## Solid Phase Extraction with Waters Sep-Pak Vac™

Most solid phase extraction (SPE) cartridges look pretty much the same. But over time, the results they produce can vary significantly. That's why scientists look for batch-to-batch and lot-to-lot reproducibility... ...guaranteed when you use Sep-Pak Vacs from Waters - the pioneer of solid phase extraction technology.



### Invest in Waters quality at an attractive price

Compare our prices and see for yourself. In addition to being the most reliable, our Sep-Pak Vacs are also one of the most economical ways to prepare your samples - manually, or automatically with robotics, or batchwise with Waters Vacuum Stations.

RELATIVE RETENTION (k)			
COMPOUND	BATCH A	BATCH B	BATCH C
DHA (No)			
PROPRANOLOL	4.0	3.8	4.0
DOKEPIN :	19.8	18.8	19.6
TOLUENE	2.2	2.3	2.2
<b>ICENAPHTHENE</b>	1.5	1.6	1.6
ROPYLPARABEN	112	1.2	1.2
RUTYLPHITHALATE	4.5	4.6	4.5







Don't miss our special offer attached !