Contributors

Molecular Dynamics (Firm)

Publication/Creation

[Chesham] : Molecular Dynamics, 1996.

Persistent URL

https://wellcomecollection.org/works/z86jz32x

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection 183 Euston Road London NW1 2BE UK T +44 (0)20 7611 8722 E library@wellcomecollection.org https://wellcomecollection.org



Not just for big budgets anymore.

Fast, quantitative autoradiography without film

You no longer need a big budget to escape the time-consuming re-exposures,

inexact quantitation, and [||UJµ||U| ||||Uyu| ; JI darkroom

hassles of film autoradiography. The compact PhosphorImager SI system uses sensitive Kodak[™] storage phosphor screens for fast, quantitative imaging of 14C, 3H, 125I, 131I, 32P, 33P, 35S, and other radioisotopes, and delivers publicationquality images in one-tenth the exposure time of conventional X-ray film. And now, new lower pricing on the PhosphorImager SI makes these advantages available to everyone.

Field-proven performance

The PhosphorImager SI system is a real Molecular Dynamics PhosphorImager. PhosphorImager systems from Molecular Dynamics are the most popular-and widely cited-filmless autoradiography systems in the world.





reply card. Or call us at (44) 1494 - 793377

There's more PhosphorImager info on our web site too-www.mdyn.com.



Sak C



ImageQuant® software (included) for Macintosh® or Windows® workstations simplifies analysis and lets you print out complete, organized records of your data and results.



IBRS PHQ D 1874 HP



Molecular Dynamics Ltd 5 Beech House Chiltern Court Asheridge Road CHESHAM Bucks GREAT BRITAIN HP5 2YZ





After exposure, slide the screen into the compact PhosphorImager SI. Scanning takes just a few minutes.

Kodak storage phosphor screens are simple to use



The Molecular Dynamics

PhosphorImager[™] SI

is now compatible with

ordinary budgets.

Please send me the PhosphorImager SI brochure.

	Molecular Dynamics							
Innovative Instrument Systems for the Life Sciences				Name		Title		
Molecular Dynamics Ltd. 5 Beech House, Chiltern Ct.								
				ny	Institution			
Asheridge Road								
Chesham, Bucks., HP5 2PX				Department				
ENGLAND								
				Address				
For fai	ster respon	nse pleas	se					
call (01494) 793377 or					Posta	Postal Code		
	1494) 793	53/1 Or						
fax to (01494) 793222 in the U.K. or call one of the numbers below				Country				
					Best time to call			
							DP3	
Argentina Australia	(01) 382-4020	Germany Greece	(02151) 83 87 0	Kuwait Mexico	43 32 555	Spain Sweden	(93) 456-06-07	
Belgium	(03) 9810-9572 (02) 2120-733	Hong Kong	01 68 95 260 (02) 578-5839	The Netherlan	(05) 341-7764	Switzerland	(08) 621-34-00 (061) 261-66-65	
Brazil	(011) 530-7833	Hungary &	(02) 370-3033	The wetherian	(34) 65 50 55 6	Taiwan	(02) 836-5831	
China	86 2 345 1100	Romania	(01) 25 10 344	New Zealand	(09) 377-0392	Thailand	(02) 424-0942	
Czech Rep.	02-20 61 01 51-2	India	(044) 434-0174	Norway	(022) 90-00-00		or (02) 435-1695	
Denmark	(043) 86-87-88	Israel	(03) 95 12 020	Poland	(022) 39-2301	Turkey	(0312) 231-3829	
Eire	(012) 95 25 01	Italy	(02) 89 12 18 08	Portugal	(01) 35 24 984		or (0312) 231-966	
Finland	(0) 804-551	Marra	or (06) 90 67 01 (02) 644-1994	Singapore	779-1919	U.S.A.	1 (800) 333-5703	
France	(1) 60 86 65 13	Korea	(02) 044-1994	South Africa	(011) 79 23 300	1		