

**Designed for comfort and use... : whatever your application : Finnpiquette  
digitals from Labsystems.**

**Contributors**

Labsystems (Firm)

**Publication/Creation**

Helsinki : Labsystems, 1991.

**Persistent URL**

<https://wellcomecollection.org/works/rk6jjjk8>


**License and attribution**


Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome  
collection**

Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

Designed for comfort and use....  
Whatever your application.



Finnpipette® Digitals from  Labsystems.


# Finnpipette® Digitals have all the features you need



**Ease of use and comfort** — Simply turn the thumb button. It clicks to the desired volume. No scales to align. Just check it on the clear, easy-to-read digital display and start pipetting. Without straining your hand. . . even after hours of use — thanks to the lightness and comfort common to all Finnpipette Digitals. . .



**Safety and durability** — The streamlined unobtrusive tip ejector ensures safety even with biohazardous materials. And also prevents accidental detipping. Durable construction and space-age plastics mean longer bench life.




**Accuracy plus in-lab recalibration** — Each Finnpipette arrives ready for use with its own accuracy validation certificate. What's more, recalibration can be carried out in your laboratory, quickly and easily.

Whatever the application you're sure to find a suitable Finnpipette Digital from our range listed below. For best results we also recommend using Finntips. Made from clear contamination-free plastic, they are available vertically packed for added convenience in boxes, or bulk packed in economically sized plastic bags.

Ordering code	volume range	scale division	accuracy %	precision %	Standard tip ordering code	package size
4027000	0.5–10 ul	0.1 ul	±5 to 1	5 to 1	Finntip 60 9400070	box of 500
4027010	5–40 ul	0.5 ul	±2.5 to 0.5	2.2 to 0.5	Finntip 60 9400070	box of 500
4027020	40–200 ul	1.0 ul	±1.1 to 0.4	0.8 to 0.4	Finntip 60 9400070	box of 500
4027030	200–1000 ul	5.0 ul	±0.9 to 0.3	0.6 to 0.3	Finntip 61 9401070	box of 200
4027040	1–5 ml	0.05 ml	±1 to 0.6	0.4 to 0.2	Finntip 62 9402050	box of 75
4027050	2–10 ml	0.1 ml	±0.3 to 0.4	0.1 to 0.3	Finntip 63 9402150	box of 40
					Finntip Universal*9400210	tray of 96

\*Ask for details of other pipettes — both single and multi-channel — in our range.

\*Fits pipettes in the 0.5–200 ul range.

 Labsystems

Manufactured by Labsystems, P.O.Box 8, SF 00880 Helsinki, Finland, Tel. +358-0-75 821, Telefax +358-0-7591 479