

**Pure DNA : improve your DNA purification with centrifugal devices from Millipore : Amicon bioseparations / Millipore.**

**Contributors**

Millipore Corporation.

**Publication/Creation**

[Place of publication not identified] : Millipore Corporation, 1999.

**Persistent URL**

<https://wellcomecollection.org/works/xh4dmn6j>

**License and attribution**

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

MILLIPORE


## More improvements...

... for your DNA applications using the unique  
MultiScreen® 96-well Separation System

- Plasmid preparation
- PCR clean-up
- Dye terminator removal

Please tick the appropriate box on the enclosed mail form.  
Full information is waiting for you!

For more information on MultiScreen  
[www.millipore.com/multiscreen](http://www.millipore.com/multiscreen)



Visit Millipore on the Internet  
[www.millipore.com](http://www.millipore.com)

Millipore, Amicon, Microcon, MultiScreen, and Ultraclean are registered trademarks  
of Millipore Corporation or an affiliated company.  
Copyright © 1999, Millipore Corporation or an affiliated company. All rights reserved.  
Cat. No. DM221EN00. Printed in France, 9/99. 99-167


**amicon**  
BIOSERPARATIONS

MILLIPORE

Special offers  
enclosed

## pure DNA


Improve your  
DNA purification  
with centrifugal devices  
from Millipore



**amicon**  
BIOSERPARATIONS

## PCR clean-up

No binding! No eluting! Just spin!



Microcon®-PCR Filter Units purify PCR products without ethanol or salts in a standard microfuge — with excellent recoveries.

**Fast** — 15-minute filtration procedure


**High recovery** — Suitable for sequencing and cloning

**Robust technique** — Separation is not affected by buffer composition, pH, or experience

**Flexible** — One device suitable for concentration and contaminant removal (primers, dNTPs, and salts)

Special offers for  
Microcon-PCR Filter Units  
and Ultrafree-DA Units.  
See enclosed mail form.

For more information  
[www.millipore.com/puredna](http://www.millipore.com/puredna)



### Simple Protocol:

1. Add PCR sample to device
2. Spin for 15 minutes at 1,000 xg
3. Invert and spin to recover sample

## DNA extraction from agarose

Fully functional DNA in one single spin!


Ultrafree®-DA Centrifugal Filter Units fragment and compress the gel, and extract DNA — all in one operation. Extruded DNA is ready for sequencing.

**Fast** — 10-minute procedure

**Highest recovery** — Recover fully functional DNA with minimal loss

**Robust technique** — Not affected by buffer, pH, or experience

**Convenient** — Handling multiple samples; no pipetting steps



### Simple Protocol:

1. Place gel slice into device
2. Spin for 10 minutes at 5,000 xg
3. Recover DNA for sequencing