

Special offer : introducing the new Millipore Pyrogard D Ultrafilter / Glenn Young.

Contributors

Young, Glenn.
Millipore (UK)

Publication/Creation

Watford : Millipore, [2002]

Persistent URL

<https://wellcomecollection.org/works/tnrfxjvd>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

Special Offer

MILLIPORE



Introducing the new Millipore Pyrogard™ D Ultrafilter

Ribonucleases (RNases) can dramatically affect results in any molecular biology experiments involving RNA and/or DNA. RNases are ubiquitous, found in human skin, body fluids and micro-organisms in air and water for example, and picogram quantities are sufficient to degrade RNA. The presence of RNases has been shown to have a negative impact on the following critical applications:

- RNA purification or isolation from cell lysate
- oligonucleotide synthesis
- RNA hybridisation
- RNA amplification such as PCR (polymerase chain reaction) and RT-PCR (reverse transcriptase-PCR)
- plasmid preparation



Pyrogard D from Millipore is a disposable ultrafilter that can be used to remove nucleases and pyrogens at the point-of-use of Milli-Q, Synergy and Simplicity ultrapure water systems.

The Pyrogard D cartridge offers a simple, cost-effective and convenient method for producing nuclease-, pyrogen-, and bacteria-free water and a safe alternative to DEPC treatment.

If your ultrapure water system is not designed to produce RNase- and pyrogen-free water, if you currently purchase bottled water, or use DEPC then try Pyrogard D by taking advantage of our 50 % discount trial offer. Just fill in the fax-back form, or alternatively post it to the address shown overleaf.

Glenn Young
Area Division Manager
Millipore (U.K.) Ltd.

P.S.: The Pyrogard D is packaged in a box of two cartridges costing £82.00. This trial offer entitles you to a 50 % discount on your first box.

Special Offer



MILLIPORE

Mr. Stephen Lowther
Wellcome Institute
Library
183 Euston Road
London NW1 2BE
United Kingdom

Fax Back Form to: (01923) 816128

Please tick the appropriate box:

- ☐ I would like to find out more about Pyrogard D

R&D Notebook RD006: Qualification of a point-of-use ultrafiltration device for the production of pyrogen-free ultrapure water

R&D Notebook RD007: Trace ion chromatography

R&D Notebook RD008: Comparison study of ultrafiltration vs. DEPC treatment for RNase removal

Pyrogard D Product Information Data sheet

- ☐ I would like more information on Millipore water purification systems
Milli-Q, Simplicity and Synergy systems literature

- ☐ I would like to take advantage of this special trial offer

Your regional Millipore Sales Representative will contact you soon to answer any additional questions on how Pyrogard D can help improve your daily work at this very affordable price:

E-mail: _____

Telephone: _____ Fax: _____

Please fax or post this form to: Anne Brown
Millipore (UK) Ltd
Units 3 & 5 The Courtyards
Hatters Lane, Watford
Herts.
WD18 8YH.
Fax: (01923) 816128

