

Introducing the revolutionary new FluorImager 575 from Molecular Dynamics : fluorescence imaging : moves you into the fast lane / Molecular Dynamics Limited.

Contributors

Molecular Dynamics (Firm)

Publication/Creation

Kemsing : Molecular Dynamics Limited, [1993]

Persistent URL

<https://wellcomecollection.org/works/vmdf6eyw>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).



Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

INTRODUCING THE REVOLUTIONARY
NEW FLUORIMAGER™ 575 FROM MOLECULAR DYNAMICS.

FLUORESCENCE IMAGING



MOVES INTO THE FAST LANE

There's a new detection technology that makes fluorescence analysis as accurate as autoradiography. And a lot faster than chemiluminescence.

The Molecular Dynamics FluorImager 575.

The FluorImager 575 directly scans samples up to 20x25 cm in less than 5 minutes. That's because there are no drying, exposure or film development steps – or radioactive materials – to slow you down.

In fact, with dual-wave-length, two-color scanning, the FluorImager 575 can run

samples and standards in the same lane. So you can get twice as much done in a fraction of the time it used to take.

Even more important, you get better results.

From gels, blots, TLC and microtitre plates, the FluorImager 575 delivers 100 times the sensitivity of traditional

fluorescence methods, and up to 4 orders of magnitude linear dynamic range. That means you can quantitate faint and dark bands without re-running samples.

Specially developed FluorKit™ reagents and standard protocols

ensure reproducibility.


And Image-QuaNT™, Molecular Dynamics'

proven analysis software, simplifies quantitation and makes it easy to transfer data to your favorite Windows and Macintosh programs.

So, if you're looking for a faster way to get quantitative results, move to the FluorImager 575. And leave the rest behind.

...

Call today for a copy of our new FluorImager 575 brochure. In the U.S., call 1-800-788-0634. In Europe, call 44 732 762565. In Australia, call (03) 810-9572. In Japan, call (03) 3976-9692.



Call us now for a no-obligation demonstration or use this leaflet as a return mailer (see reverse)



Molecular Dynamics Limited, 4 Chaucer Business Park, Kemsing, Sevenoaks, Kent, TN15 6PL, UK. Tel: 0732 762565 Fax: 0732 763422

©1993, Molecular Dynamics. FluorImager™, FluorKit™, and ImageQuaNT™ are trademarks of Molecular Dynamics. MD003

FOLD 3 AND TUCK

PLEASE
PLACE
STAMP
HERE

FOLD 1

TO: MOLECULAR DYNAMICS LTD.,
4 Chaucer Business Park,
Kemsing,
Sevenoaks,
Kent,
TN15 6PL,
UK.

FOLD 2

- ☐ Please send me detailed information on the FluorImager
- ☐ Please have a representative contact me for a demonstration
- ☐ Please send information on the Molecular Dynamics Product Range

NAME: _____ TELEPHONE: _____

DEPARTMENT: _____ FAX NO: _____

INSTITUTE: _____

ADDRESS: _____

TOWN: _____

POSTCODE: _____

