

**"Spiromol" solutions.**

**Contributors**

Allen & Hanburys

**Publication/Creation**

London : Allen & Hanburys, [1925]

**Persistent URL**

<https://wellcomecollection.org/works/k6jgky4z>

**License and attribution**

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome  
collection**

Wellcome Collection  
183 Euston Road  
London NW1 2BE UK  
T +44 (0)20 7611 8722  
E [library@wellcomecollection.org](mailto:library@wellcomecollection.org)  
<https://wellcomecollection.org>

# "Spiromol" Solutions

Trade Mark

*For Oro-Nasal Inhalation in Phthisis,  
Bronchitis, Nasal Catarrh, etc.*

Attention has recently been directed to the great benefit that may be derived in Phthisis from inhalation of such antiseptics as Iodine, Carbolic Acid, Formaldehyde and Creosote, especially when combined, where cough is troublesome, with Spirit of Chloroform or Ether.

## *Composition :*

- "Spiromol" No. 1  
Ethereal Tincture of Iodine, Carbolic Acid, Creosote  
(Dr. Coghill's formula).
- "Spiromol" No. 2  
Ethereal Tincture of Iodine, Carbolic Acid, Thymol.
- "Spiromol" No. 3  
Carbolic Acid, Creosote, Iodine, Ether, Chloroform  
(Dr. Lees' formula).
- "Spiromol" No. 4  
Formaldehyde, Chloroform, Menthol, Oil of Pine  
(Dr. Muthu's formula).
- "Spiromol" No. 5  
Formaldehyde, Ether, Oil of Pine (Dr. Wilkinson's  
formula).
- "Spiromol" No. 6  
Iodine, Creosote, Eucalyptus, Terebene, Chloroform,  
Menthol, Ether (Dr. Macdonald's formula).

## *Instructions for Use :*

Place in the inhaler sufficient of the Solution to cover the end of the elongated tube, then adjust the stopper so that the air enters freely, and the inhaler is ready for use, the black dots being in line. Place in nostril and draw long breaths causing air to bubble through bottle. Opposite nostril can be closed with finger if wished. Breathe out through mouth. Put in some drops of fresh inhalant frequently.

## *Price of Solutions :*

1 oz. 2s. 6d.      2 oz. 4s. 6d.

# All Glass Inhaler

Patent No.

174,540

For the  
Prevention and  
Cure of  
Influenza, etc.



Price :  
In strong card carton,  
2s. 6d.



SCALE  $\frac{1}{2}$

Designed by  
Dr. R. C. Macdonald,  
INVERNESS.



The "Macdonald" Inhaler may be sterilised by boiling, and by a quarter turn of the stopper it is securely sealed, and may be carried in pocket or hand bag without risk of leakage.

This is a new pattern inhaler so constructed that by nasal or oral inspiration air is drawn into the main body of the bottle-like apparatus through a small opening in the neck. It then passes down the central tube, through the inhalant, thus becoming saturated with the medicaments. The medicated air then enters the nozzle and is inhaled by the patient through the nose or mouth as required. The inhaler gives a highly saturated vapour with drugs of a volatile nature. It has been found beneficial in cases of Phthisis, Laryngeal Tuberculosis, Nasal Catarrh, Rhinitis, Suppuration of accessory sinuses, antrum or frontal, Naso-Pharyngeal Catarrh, Asthma and Bronchial affections, and as a preventive of Influenza.

For particulars of Solutions for use with Inhaler see over.

## Allen & Hanburys Ltd.

Manufacturers of Surgical Instruments & Appliances,  
48 Wigmore Street, LONDON, W.1.

Telephone 6476 Mayfair (four lines).

Telegrams: "Orthopedic, Wesdo, London."