

For preserving eggs : silicate of soda or water glass / Peter Leech & Son.

Contributors

Peter Leech & Son.

Publication/Creation

Whitehaven : Peter Leech & Son, [between 1910 and 1919?]

Persistent URL

<https://wellcomecollection.org/works/p9vmx6uk>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome
collection**

Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

**DIRECTIONS FOR
PRESERVING EGGS.**

TO make the Preserving Liquid, add to the contents of this tin from two to three gallons of boiling water and stir until dissolved, then set aside until quite cold.

Place the Eggs carefully in a wooden or earthenware vessel and pour the Preserving Liquid over them until well covered. Tie a cloth over the vessel to keep out the dust and dirt. The vessel to be kept in a cool place.

New Laid Eggs treated by this method will maintain their milky freshness for over six months.

FOR PRESERVING EGGS.

Silicate of Soda
—OR—
Water Glass.

Peter Leech & Son,



WHITEHAVEN.