

Pifco vibratory massager : the way to glowing health... soothes away pain and fatigue / Pifco.

Contributors

PIFCO Ltd.

Publication/Creation

[Manchester] : Pifco, [1948?]

Persistent URL

<https://wellcomecollection.org/works/b9w65ud7>

License and attribution


Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome
collection**

Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

PIFCO *Vibratory Massager*

*The way to
glowing
health*



The PIFCO Massager enables complete massage to be enjoyed economically in the privacy and convenience of the home, with all the benefits of an expensive massage.

For Rheumatic ailments, stiffness and other aches and pains the PIFCO Massager is invaluable. For toning up the body and for slimming.

A daily face massage promotes fresher and lovelier skin, soothes away lines of tiredness, brings inner loveliness to light.

Massage helps rid the scalp of Dandruff, brings new life to tired hair follicles.

Give your family and yourself that zest for life that goes hand in hand with a glowing healthy body.

The PIFCO Massager is beautifully styled in ivory plastic with intensity regulator and supplied complete with four applicators and comprehensive instruction manual. 6 ft. circular black and white braided flex and adaptor fitted.

Soothes away pain and fatigue

Sponge Rubber Applicator



Very refreshing
and stimu-
lating.
Helps
keep a
clear and
radiant skin.

Spiked Rubber Applicator



Stimulates the
circulation
round the hair
roots. Helps
keep the hair
healthy and
scintillating.

Flat Rubber Applicator



Alleviates Rheu-
matic and
kindred aches
and pains.
Helps disperse
fatty accu-
mulations.

Hard Applicator



Athletes use the
Pifco Massager
to tone
up muscles.
Helps
avoid
stiffness.



The PIFCO Massager
complete with four Appli-
cators and comprehensive
instruction manual.
No. 1056. For 200/250
Volts. A.C. only.

85/6 complete

Also made in 100/120 volts.

Vibratory Massager Carrying Case

Velveteen lined, with shaped platform to carry
the Massager and four Applicators. Exterior
rexine covered in red and black, with hinged lid
and leather carrying handle.

Dimensions: 10 $\frac{3}{4}$ in. \times 8 in. \times 3 $\frac{1}{4}$ in.

Carrying Case only. No. 1057 47/11 complete

