

Polybactrin bacitracin-neomycin-polymyxin : antibiotic powder spray.

Contributors

Calmic Ltd.

Publication/Creation

[between 1960 and 1969?]

Persistent URL

<https://wellcomecollection.org/works/sv36r6zm>

License and attribution

Conditions of use: it is possible this item is protected by copyright and/or related rights. You are free to use this item in any way that is permitted by the copyright and related rights legislation that applies to your use. For other uses you need to obtain permission from the rights-holder(s).

**wellcome
collection**

Wellcome Collection
183 Euston Road
London NW1 2BE UK
T +44 (0)20 7611 8722
E library@wellcomecollection.org
<https://wellcomecollection.org>

polybactrin bacitracin-neomycin-polymyxin


(T.S.A. LICENCES 2/308/0126)

ANTIBIOTIC POWDER SPRAY
CALMIC LIMITED. CREWE. ENGLAND.

polybactrin bacitracin neomycin polymyxin

(T.S.A. LICENCES 2/308/0126)

ANTIBIOTIC POWDER SPRAY
CALMIC LIMITED. CREWE. ENGLAND




SHAKE FREQUENTLY
polybactrin
bacitracin
neomycin
polymyxin
POWDER SPRAY
T.S.A. LICENCES 2/308/0126
STERILE

INSTRUCTIONS

1. Use and during use.
2. Maximum storage period.
3. Do not puncture or use press applicator.
4. Do not autoclave.
5. Do not puncture or use press applicator.
6. Do not autoclave.

Manufactured in England by
CALMIC LIMITED · CREWE

USE THIS WAY UP



SHAKE FREQUENTLY
polybactrin
bacitracin
neomycin
polymyxin
POWDER SPRAY
T.S.A. LICENCES 2/308/0126
STERILE

INSTRUCTIONS

1. Use and during use.
2. Maximum storage period.
3. Do not puncture or use press applicator.
4. Do not autoclave.
5. Do not puncture or use press applicator.
6. Do not autoclave.


Manufactured in England by
CALMIC LIMITED · CREWE



Each batch of Polybactrin passes through 7 laboratory checks to ensure sterility and potency before being released for distribution.

BROAD SPECTRUM ACTIVITY

STERILE



FORMULA Each canister contains:
Neomycin Sulphate ... = 495 mg. base
Polymyxin B Sulphate ... 150,000 units
Zinc Bacitracin ... 37,500 units
Pressurised with dichlorotetrafluoroethane
70% and dichlorodifluoromethane 30%.
(109 g. approx.)
Maximum storage temperature 49°C (120°F).

POLYBACTRIN

The effectiveness of Polybactrin against a wide range of pathogens has been effectively proved in clinical studies extending over a period of seven years. In using Polybactrin,

Surgeons may be sure that they have a safe and potent topical application of antibiotics for use as routine prophylaxis and treatment in all surgical procedures.

CALMIC LIMITED, Crewe, Cheshire Tel: 3251 (7 lines)